

SIEMENS MAGNETOM Prisma syngo MR D13D

\\USER\NTU\Brain\NTU_CCY\ep2d_bold_p2_64_3x3x3_186
 TA:9:29 PAT:2 Voxel size:3.0x3.0x3.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	45
Dist. factor	0 %
Position	R2.0 P2.3 H5.5 mm
Orientation	T > C-6.6 > S3.2
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	3000 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HE1-4

Contrast

MTC	Off
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	186
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	On
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	45
Dist. factor	0 %
Position	R2.0 P2.3 H5.5 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P

System

Body	Off
HE1	On
HE3	On
NE1	Off
HE2	On
HE4	On
NE2	Off
SP5	Off
SP6	Off
SP7	Off
SP8	Off
SP1	Off
SP2	Off
SP3	Off
SP4	Off
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	R2.0 P2.3 H5.5 mm
Rotation	1.40 deg
R >> L	192 mm
A >> P	192 mm
F >> H	135 mm
Frequency 1H	123.255624 MHz
Correction factor	1
SincRFPulse 1H	270.091 V

Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
-----------------	------

Inline

Distortion correction	Off
-----------------------	-----

Sequence

Introduction	On
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2442 Hz/Px
Free echo spacing	On
Echo spacing	0.49 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HE1-4
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	On
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	6
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off